

### MEDICAL STANDARDS STABLE SUPPORT

Medical molecular sieve oxygen producer 8F-5AW [5L machine]



### SIX FUNCTIONS, NEW UPGRADE



8F-5AW

90%-96%

Medical oxygen concentration



F 4



Intelligent remote control

A|B

A/B Oxygen production and atomization

15.5<sub>KG</sub>

Small in size and light in weight




Streamline design is more

fashionable



# DOUBLE LAYER FILTRATION OXYGEN MORE PURE AND CLEAN



### EPL COMPRESSOR CONTINUE TO PROVIDE YOU WITH GREAT MOTIVATION

Strong wind pressure makes the heat dissipation performance of the machine significantly improved, Effectively extend the service life of oxygen generator.



### LOWER NOISE AND MORE COMFORTABLE

Admission air silencer, exhaustsilencer, Stereo cycle noise reduction technology to give you a meditation world



### PROFESSIONAL ATOMIZATION CAN MATCH VENTILATOR

One machine dual use Oxygen production atomization switch at will



Oxygen atomization Oxygen can be inhaled is used at the same timeConstant oxygen same time Constant concentration

by two persons at the oxygen concentration

### REMOTE CONTROL FUNCTION IS MORE CONVENIENT

The wireless remote control can be adjusted as it pleases





# LIGHT AND PORTABLE EASY TO MOVE

NEW 5L OXYGEN GENERATOR, WEIGHT IS CHANGED.

Only 15.5 KG


### REMOVABLE FILTER FELT



## ONE-KEY START MORE CONVENIENT TO USE



## FAULT DETECTION FUNCTION



UUWell 鱼跃

- Current overload or loose connection
  - shut down
- Compressor high temperature
  - shut down
- Pressure, cycle fault shut down
  - -ṁ-alarm ⊝ shut down
- Low oxygen
   concentration
  - ப் alarm



#### Product Specification

Product model: 8F-5AW

Power:500VA

The oxygen flow rate:1-5L/min

After-sales:The machine is guaranteed for 1 year,the compressor for 3 years

Main raw material:UOP molecular sieve Noise:Standard values is 60dB(A),Yuwell measurement value is 48dB(A)

Net weight: 15.5KG

Oxygen

concentration:0.5-5L/min,93%±3%

Dimensions:39\*24.5\*50cm

Compressor:Epl

#### Matters need attention

#### Special warning

To prevent power failure or possible failure of the oxygen generator, Those in urgent need of oxygen and those in critical condition must be equipped with alternative feeding devices (such as: Oxygen cylinder, oxygen bag, etc.)

The equipment is suitable for oxygen supplement, not as life support and life extension supplies.